


USER MANUAL

Athletic Pro Massage Chair




Please read all the instructions before using the chair.

1. Unpacking the chair

The chair is delivered in a single box. Get at least two people to carefully unpack the main body and accessories from the carton box. Ensure all the accessories are present and accounted for:

- ✓ Leg rest
- ✓ Main body
- ✓ Front cover & pillow
- ✓ Remote control
- ✓ Power cable



Dimensions

Packing size (L×W×H): 46" x 31" x 50"

Product dimensions (upright): 52" x 28" x 46"

Product dimensions (reclined): 67" x 28" x 38"

Net Weight: 185 lbs


Gross Weight: 222 lbs

Door clearance: 30"


2. Installing the massage chair

Step 1. Pad and leg rest



- Put the leg rest assembly in front of the seat near the seat front part.
- Put the pad on top and connect the air hose to the main body of the chair.
- Zip the pad on each side and on the top. Zip the pillow on top of the pad.
- Connect the air hose power wires from the leg rest to the main body.
- Make sure you do not wrap the wires with the air hose or block the air hose.
- Drop in the leg rest into the 2 metal latches / hangers.
- Make sure the level limiting pin goes over the hanger to block the leg rest when pushed up.
- Uncover the velcro and paste it to the seat, under the cushion.
- Paste the second velcro, on the back of the cushion to the top of the leg rest.



Air hose and
electrical wiring




Latches/hangers and the level limiting pin positions



Step 2. Install the remote control

- Retrieve the remote controller pocket. It is folded in the interior of the armrest.
- Put the remote controller in the pocket.
- Make sure the wire goes under the airbag in the armrest.



Remote control installation

Step 3. Positioning the massage chair

Choose your ideal installation location by taking into consideration the following requirements.

⚠ ATTENTION

- DO NOT install the chair in a high moisture environment or on an unlevelled floor.
- DO NOT position the chair next to a heat source, open flame or in direct sunlight.
- DO NOT place or run the massage chair over the power cord.

Keep the massage 3ft away from TV, radio or other equipment to avoid potential signal interference.

Allow for clearance at least 2" from the back to the wall and 20" in front of the leg rest.

Step 4. Moving the chair around

The chair has 2 wheels on the back and can be moved around by one or two persons. To safely move the chair push down on the backrest and tilt it 30°.


⚠ CAUTION

- Make sure the power cord is not connected while moving.
- Hold the backrest firmly, or you may drop the chair.
- Do not move the chair while a user sits in it.
- Soft floors may be scratched by the wheels. Place a soft mat or carpet underneath.

Step 5. Plug in the chair

Finish reading the manual before making this step.

- Connect the power cord to the power box in the back of the massage chair.
- Turn ON the power switch next to the power cable.
- To start the chair press the ON button on the remote control or the ON quick control button.
- Turn the power OFF and disconnect the power cable when you know the massage chair will not be used for extended periods of time.



3. Technical Specifications

Brand:	Rebalanse
Model:	Athletic Pro
Input voltage:	110-120V, 60 Hz
Power:	180 W
Airbags:	26
Rollers:	Yes
Massage mechanism:	2D
Manual massages:	Yes
USB port charging	Yes
Foot massage roller:	Yes
Foot massage airbag:	Yes
Armrest/Hands massage:	Yes
Zero gravity:	Yes
Full body heating:	Yes
Back heating:	Yes
Cover material:	PVC Leather
FDA Classification:	Class I Medical Device
Track Type:	S & L track
Track Length:	53"
Colors:	Navy Brown, Cream White
Bluetooth:	Yes
Body scanning technology:	Yes
Speakers:	Yes
Session preset duration:	10, 20, 30 min
Wheels:	Yes, 2
Remote type:	LCD remote + quick control buttons
Packing size (L×W×H):	46" x 31" x 50"
Net Weight:	185 lbs
Gross Weight:	222 lbs
Upright dimensions:	52" x 28" x 46"
Reclined dimensions:	67" x 28" x 38"
Door clearance:	30"

4. IMPORTANT SAFETY INSTRUCTIONS

Please read this manual carefully to ensure the correct installation and operation of the massage chair. The Safety Instructions are provided as both Warnings and Cautions for important safety information. Pay close attention to both Warnings and Cautions as they indicate essential safety precautions.

⚠️ WARNINGS

Warnings indicate potential risks of severe injury or harm to the body.

- ⚠️ A. People with conditions such as osteoporosis, vertebral fractures, acute painful diseases (sprains, muscle strains), heart disease (including pacemaker users), cancer patients, perceptual disorders due to circulatory disturbances (e.g., diabetes), abnormal spine, crooked spine, protruded intervertebral disc, pregnant women, those who have just completed childbirth, and individuals with poor health should not use the massage chair.
- ⚠️ B. Individuals with thrombosis, acute varices, various skin inflammations, unable to feel temperature levels, high fever (acute inflammation), and generally individuals over 80 years of age should avoid using the massage chair as it may cause accidents or discomfort.
- ⚠️ C. People with impaired consciousness, mental disabilities, or those unable to operate the chair independently should refrain from using it. It is also prohibited for infants, individuals who cannot take care of themselves, and those under the influence of alcohol.
- ⚠️ D. Do not use the massage chair when you are wet. It should not be used by individuals who are not sensitive to heat or who have allergies to heat.
- ⚠️ E. Individuals experiencing any of the above-mentioned symptoms should avoid using the massage chair or consult a doctor before usage. Users should exercise discretion when using the chair.
- ⚠️ F. Prolonged usage (more than 1 hour) is not recommend, nor is sleeping in the chair!

⚠️ FIRE AND ELECTRIC SHOCK PRECAUTIONS

- ⚠️ Avoid running the chair without any load for an extended period.
- ⚠️ Never leave the chair unattended while it is plugged in and working.
- ⚠️ Ensure that the power cord, plug, and socket are not broken or loose.
- ⚠️ Use only the rated voltage and avoid using other electric power sources.
- ⚠️ Connect the chair's cord, which includes a grounding conductor, only to grounded sockets.

Note: If the fuse is burnt out, turn off the switch, unplug the power, and open the lid to access the fuse in the switch box. Replace it with the same-size fuse and close the lid. Alternatively, consult an electrician for fuse replacement.

5. ⚠ OPERATING PRECAUTIONS ⚠

- DO NOT attempt to move the chair while it is in operation.
- DO NOT place your hands or feet in the gaps between mechanical parts.
- DO NOT get into the chair when the footrest is lifted.
- DO NOT use the chair directly on bare skin as it may cause skin irritation.
- DO NOT use electric blankets or other heating products on the chair.
- DO NOT stand on the chair while it is in use.
- Make sure not to drop or insert any objects into the massage chair.
- Avoid storing hard objects in your pocket while getting a massage.
- Keep pets and babies away from the chair, especially the area under the leg rest.
- Please remove any hair ornaments before using the massage chair.
- Avoid unplugging the chair or turning off the power during operation.
- The chair is intended for indoor use only and should not be used outdoors.
- Place the chair away from heat sources, water, dust, corrosive, inflammable environments.

6. Operation Guide

BLUETOOTH PAIRING

By default, when the massage chair is ON, the Bluetooth is on, accepting connections. No password is required for pairing. The chair speakers will be used to play sounds from the connecting device. Their maximum volume level is 80 decibels (dB), roughly equivalent to a louder conversation. Bluetooth can be deactivated from the ⚙ Settings menu.

REMOTE CONTROL / QUICK CONTROL

Buttons on the remote:

On/Off: press this button to turn on/off the massage chair.

Auto programs: press this button to enter the automatic massage programs. By default they run for 20 minutes, and can be adjusted to 10 or 30 minutes. The automatic programs are:

Recover	Relief	Stretch	Music Sync	
Renew	Calming	Special	Refresh	Energy

Airbag program: press this key to start the full-body, upper body, or lower body air pressure massage functions in sequence.

Emergency: press this button to immediately stop all massage functions. Afterwards, only the shutdown button will be accessible for selection.

Menu: press this button to access the massage chair's options.

Auto	Manual	Airbag	Roller
Heating	Timing	Setting	Stop

The menu key needs to be used in tandem with the navigation keys.

Method: press this key to switch the massage method/technique: kneading, tapping, knocking, kneading & tapping, and shiatsu. Press OK to confirm the selection of menu items.

Right, Left, Up, Down: are the navigation keys used within the menu. In the menu mode, you can use the up and down direction keys to select menu items and the left and right direction keys to navigate between menus. Additionally, the up and down direction keys help you adjust the position of the massage mechanism during shoulder height detection and manual massage mode.

Heating: press this key to turn on and off the waist heating function. When turned on, the corresponding heating symbol turns red, and when turned off, the corresponding heating symbol disappears.

Position: when you are in the manual massage mode, you can press this button to select Whole, Partial and Point.

Width: when you are in the manual massage mode and using the tapping, knocking or shiatsu massage techniques, you can press this button to adjust the massage width: narrow, middle, and wide.

Intensity: use this key to adjust the intensity of the air pressure on 5 levels.

Speed: when in the manual massage mode, press this key to adjust the 5 massage speeds.


Foot roller: press this key, and the sole roller will start the massage. You can choose Quick, Soft or to Stop.

Leg rest up/down: By pressing and holding these buttons, the leg rest will start to rise/descend.

Back up & Leg up/down: By pressing and holding these buttons, the backrest and leg rest will ascend/descend. Releasing the button will halt the movement of both the backrest and leg rest.

Zero gravity: press this button to bring the chair to the zero gravity position. Then press this button again to reset the chair to the initial working state. The zero gravity icon on the main interface will show the current zero gravity status once the massage starts.

REMOTE CONTROL FUNCTIONS



QUICK CONTROL BUTTONS

7. Storage and Maintenance

PARTS MAINTENANCE


No maintenance is required for the components of the massage chair. It does not need oiling.

CLEANING

Taking care of your massage chair ensure optimal performance. Regular cleaning will help maintain its appearance and hygiene. Follow these simple steps to clean your massage chair:

A. Synthetic Leather and Plastic Cleaning:

- Use a tidy, soft, and dry cloth to clean gently.
- Avoid using products that contain gasoline, alcohol, or other chemicals.
- Dip a soft cloth into a 3-5% neutral detergent solution. Squeeze to remove excess liquid and gently clean the stains on the surface. Softly dry the area with a dry cloth. Allow it to air dry naturally.



B. Cleaning Cloth Material:

- Use the same neutral detergent solution and squeeze the cloth to remove excess liquid.
- Gently clean the material surface. For stubborn stains, lightly brush the dirty parts with weak neutral detergent (and avoid excessive brushing).
- Excessive cleaning may cause damage to the cloth material.
- Allow it to dry naturally. DO NOT use an air blower.

Important Safety Precautions

- Before cleaning, always unplug the massage chair from the electricity supply.
- Do not touch the plug or power socket with your wet hands.
- Ensure your hands are dry when plugging the chair back into the socket.
- Failure to follow these instructions may result in electric shock or burns.

STORAGE

Disconnect the power cord and keep the chair out of the sunlight or heat sources. Do not stack heavy items on top of it.

8. Troubleshooting

A. If the massage chair doesn't work after power is turned on:

- Make sure the control switch is turned on.
- Ensure that a massage function has been selected.
- Check if the plug and socket are properly connected. They shouldn't be in poor contact.
- If you suspect the wire or plug is damaged, get in touch with the manufacturer or maintenance personnel for a check and possible replacement.
- In case the fuse has burned out, use the replacement fuse or one of the same specification.
- If you suspect something is wrong with the inner circuits, ask a maintenance person to check and repair it.

B. Hearing a beep sound inside the massager while it's working is normal. It's the sound of the air pump, motor, and mechanical parts.

C. The different height between the left roller and the right roller is also normal. The massage rollers work on shifts and don't require any action.

D. If the sound becomes much louder while the chair is working:

- It could be due to the chair working for a long time. In this case, turn off the power and let the massager rest for an hour before you use it again.
- If the issue persists, it could be because the inner parts have been worn away due to long-term use. If this is the case, reach out to maintenance personnel for a check and possible repairs.

E. If the massage function suddenly stops while working:

- Check if there's a power outage and try connecting the electricity again.
- When the chair works for long periods of time it may overheat and activate the internal protection (stops working). Just let the massager rest at least an hour before you use it again.

F. If neither the backrest nor the leg rest cannot raise/recline:

- Check for any obstacles and remove them.
- Check to see if the connecting wires/plugs are loose.

G. If the electric wire or plug gets heated:

- This could be because the chair is overloaded. To resolve this, turn off the power and let the massager rest for one hour before the next use.
- If you notice any other abnormal case, stop using the chair and reach out to maintenance personnel for help.

CONTACT US: Visit www.Rebalanse.com for up to date contact information.